

## PR0735 Acetabulum

In an effort to assist you with your hip prosthesis courses in the future, we have developed a better consistency in the material used for the right hemi-pelvis with acetabulum malunion. The new material is not brittle and grates easily with enough elasticity. The result is a perfect setting for the hip prosthesis socket. Thanks to the integrated holder, the model can be securely mounted allowing for more precise work.



### **PR0735 Hemi-Pelvis malun.**

Description: Right Hemi-Pelvis with malunion at the Acetabulum and posterior holding block

Dimensions: Acetabular diameter: 49mm

Material: Solid foam, light dense cancellous bone

Fracture: Not fractured

Remarks: Used for hip prosthesis

